

**RNA Transcription, Purification, Characterization,
And
NMR Sample Preparation**

**.....'How To' series, 2004
January 13th**

Steps Involved

Synthesis of RNA sequence

(*in vitro* T7 Polymerase-based RNA transcription)

Purification of Synthesized RNA

(denaturing PAGE, Electro-elution)

Characterization

(PAGE, UV Melting, NMR)

NMR Sample Preparation

(Buffer choice, Lyophilization vs Speed-Vac, Annealing)

T7 RNA Polymerase-based *in vitro* Transcription

To produce large quantities of RNAs for NMR

Limitations:

Requires 5'-end to be purine-rich!
and may produce 3'end heterogeneity!

T7 bacteriophage RNA polymerase is single polypeptide chain enzyme, extremely efficient and allow RNA synthesis in large quantity

Can be obtained commercially, but becomes expensive for large scale RNA synthesis, thus has to be expressed and purified *in-house*

T7 RNA Polymerase Information


Over-expression plasmid (pAR1219) developed by Studier et al. It has gene encoding T7 RNA polymerase under control of lac UV5 Promoter. Plasmid also carries lac Iq gene encoding lac repressor. Production of T7 RNA polymerase is induced by adding IPTG.

T7 RNA polymerase is easily proteolysed by protease OmpT and E.Coli BL21 strain cells lacking this enzyme are recommended host for the plasmid.

T7 RNA Polymerase-based *in vitro* Transcription of RNA Oligonucleotide of your interest

RNA sequence

5'-pppGGAGGCUCAGUCUGGAA-3'



Transcription Cocktail

Template DNA duplex (0.01mM)

rNTPs (4mM each)

MgCl₂ (5 times NTP conc., 20-30mM range)

5X Transcription buffer

(contains Tris-HCl, pH 8.1, 20mM spermidine, Triton X-100,
100mM DTT, PEG)

T7 Polymerase Enzyme (1mg/ml)

Sterile water

TOTAL REACTION VOLUME 10-20mL

Incubation for 2-4 Hrs.

Component	[i]/mM	Vol.add/ml	[f]/mM
Top Strand	0.02	0.60	0.0006
Bottom Strand	0.02	0.60	0.0006
ATP	100.00	0.80	4.00
CTP	100.00	0.80	4.00
GTP	100.00	0.80	4.00
UTP	100.00	0.80	4.00
MgCl2	200.00	3.00	30.00
T7 RNAPol	0.25	0.20	0.0025
5x Trans. Buffer	5x	4.00	1x
H2O		8.40	
	Final Volume	20.00	
Edta	250.00	2.40	30.00
NaAcetate	3M	2.00	.10x rxn vol.
EtOH chilled		50.00	2.5x rxn vol.

TIME		8:15	9:15	10:15	11:15	0:15	13:15	14:15	Totals
	dilution uL (in 1mL)	10	10	10	10	10	10	10	
#1	Cell								
	Vol. (uL)	1100	900	940	990	980	950	960	6820
	abs.	0.03	0.06	0.034	0.049	0.018	.01	.003	!
	OD amount	3.3	5.4	3.196	4.851	1.764	0.95	0.288	19.749
#2	Cell								
	Vol. (uL)	650	670	700	780	970	940	1160	5870
	abs.	0.021	0.046	0.04	0.038	0.005	.005	.006	!
	OD amount	1.365	3.082	2.8	2.964	0.485	0.47	0.696	11.862

Effective Variables in Transcription Reaction

- 1. DNA Template final concentration 200nM-1uM**
(best results, generally, around 500nM)
- 2. rNTP concentration (2-4mM)**
(can be lowered for specific NTPs)
- 3. Magnesium Concentration** *(most important!)*
- 4. T7 RNA Polymerase concentration (0.5-2%)**
- 5. Incubation time**

Once Transcription incubation is over...

1. Quenching the reaction by EDTA 1.5 times total Mg conc.
2. PPTN. By Ethanol/3M 5.3pH NaoAc. O/N at -20C
3. Spinning at 12g for 15min
4. Wash with 70% EtOH (chilled)
5. Spin again.
6. Save the pellet. Let it dry at RT.
7. Add 6X BPP/XCC Dye in 7M Urea.
8. Heat in 7M Urea to dissolve

Load onto 20% PAGE Denaturing. O/N

600-1000V Constant (better if you pre-run for 1hr or so)

Next Day....

1. Cut the band using UV shadowing
2. Electro-elute using 1000MW cut off membrane. (0.5x TBE)
3. Check OD at 260nm of the collected sample to quantitate RNA.
4. (USE Extinction coefficient calculated from web site)
5. BioPolymer Calculator paris.chem.yale.edu/extinct.html
5. PPT again using EtOH/NaOAc

Dialysis

6. Dissolve and dialyze against
 - Water/EDTA(10mM) 24Hrs (1L volume)
 - Water/1mM EDTA (24Hrs)
 - Water only (24Hrs.)

Collect all liquid from membrane and again take OD to measure RNA quantity- Store at -20C.


Characterization:

Take 10-20uL containing ~0.05ODs and run denaturing 10-20% PAGE with known markers to check the integrity (single band) and size (MW) of your RNA.

Trying to make Hairpin...it could also make duplex!


Check with Native PAGE in TBE or TBM. Run at 4C
Also, UV-melting can give information about the melting temperature, Buffer effects, etc.

Finally, NMR (1D Imino and TOCSY) to further check the RNA Integrity (conformation and sequence)





**Characterization of SL1 RNA Dimers and Monomers
in 10% PAGE in TBE at RT**

SL1-wt RNA melting in 0.1mM Mg-containing Buffer Annealed at RT, 37°C, 50°C and 55°C



SL1-YB O/N Incubation with NCp7





TOCSY NMR of RNA showing U and C H6-H5 signals

Relevant References

- 1. Synthesis of small RNAs using T7-RNA Polymerase
Milligan and Uhlenbeck, Methods in Enzymology, Vol 180, 1988,
Page 51-62**
- 2. RNA Oligonucleotide Synthesis using T7 RNA Polymerase and
Synthetic DNA templates, NAR, Vol 15(21), 1987, 8783**
- 3. Synthesis and Purification of Large amounts of RNA oligos
Wyatt, Chastain and Puglisi, Biotechniques, Vol 11(6) 1991, 764**
- 4. NMR of Macromolecules: A practical approach Ed. GCK Roberts
Chapter 2, 7-34.**
- 5. RNA: Protein Interactions: A practical approach,
Edited by C.W.J. Smith, Chapter 2- Preparation of
RNA:Protein complexes for X-ray and NMR.Page:37-72.**